

Experience a Total Agilent Solution

Green vacuum solutions for your lab




Agilent Technologies offers more than 7,000 vacuum products, including vacuum pumps and leak detection products, vacuum indicators, components and accessories.

MAKE YOUR LAB SUSTAINABLE!

One stop supplier for all your laboratory instruments, consumables and vacuum products.

- Go Green and create a better working environment
- Lower your lab costs
- Get access to vacuum expertise
- Improve quality


















Agilent Vacuum and Leak Detection for Sustainability

IDP Scroll Pumps:

make your workday more pleasant: an innovative scroll design reduces noise and vibration.

Lower your cost of ownership:

With Agilent oil-free vacuum pumps, you won't have to worry about oil/filter replacements, or hazardous waste disposal of used oil.

Primary Vacuum Pumps for GC-MS	 IDP-3 Dry Scroll Pump	 DS 102 Rotary Vane Oil Pump	 DS 202 Rotary Vane Oil Pump	 IDP-10 Dry Scroll Pump		
Primary Vacuum Pumps for LC-MS	 DS 402 Rotary Vane Oil Pump	 MS 40+ Rotary Vane Oil Pump	 TS 800 Dry Scroll Pump VFD	 MS 120		
Primary Vacuum Pumps for ICP-MS	 MS 40S Rotary Vane Oil Pump					
Other Vacuum Solutions for Your Lab	 Turbo Pumps	 Diffusion Pumps	 Leak Detectors	 Vacuum Gauges	 Vacuum Valves	 Vacuum Components

Learn More:
www.agilent.com/chem/vacuum

Check out the right Agilent Primary Vacuum Pump for each Lab Instrument

Agilent Vacuum Pump Model	Pump Type	Part Number	Material Description	Agilent System
IDP-3	Dry Scroll	IDP3D21 (24V DC)	IDP-3 replacement pump	SQ GCMS
		X1699-64015	IDP-3 upgrade kit EU	
		X1699-64016	IDP-3 upgrade kit US	
		X1699-64017	IDP-3 upgrade kit China	
		X1699-64018	IDP-3 upgrade kit Japan	
IDP-10	Dry Scroll	X3810-64010	IDP-10 replacement pump	TQ GCMS
		X1699-64147	IDP-10 upgrade kit	
IDP-15	Dry Scroll	X3815-64010	IDP-15 replacement pump	Q-TOF GCMS
		X1699-64146	IDP-15 upgrade kit	
TriScroll 800	Dry Scroll	X3848-64800	TS 800 scroll	LC/QTOF funnel, LC/QTOF
MS 40+	Rotary Vane	9499225M004	MS 40 for LCMS	All LCMS, Ultivo
		9499225M011	MS 40 for Ultivo	
MS 40S	Rotary Vane	X3705-64010	MS 40 for ICPMS	ICPMS Series
MS 120	Rotary Vane	X3702-64011	MS 120	TQ LCMS
DS 202	Rotary Vane	9499320M012	DS 202 dual stage Rotary vane pump	Q-TOF GCMS
AX-65	Diffusion pump	X3901-80500	AX-65 W50,90V,250W,Std CC	SQ GCMS
		X3901-80501	AX-65 NW50,165V,250W,Std CC	


Accessories & Spare Parts

Part Number	Material Description
IDP-3 TS	IDP-3 Tip Seal Kit
X3807-67000	IDP-10 scroll pump tip seal kit
X3815-67000	IDP-15 scroll pump tip seal kit
PTSS0800TS	TS 800 scroll pump tip seal kit
X3848-64850	TS 800 scroll pump module for 6550
9499390	DS pumps oil bottle AVF Silver
X3760-64005	MS 40 pumps oil bottle AVF Gold
X3760-64007	MS 40 pumps oil bottle AVS Platinum

Part Number	Material Description
9499342M003	Exhaust replacement cartridge
9499375	DS 3/8 oil drain kit
9499392M006	DS exhaust filter NW25
9499396M003	DS power cable
X3901-24000	AX-65, 250W, 90V Heater harness W Z7523
X3901-24001	AX-65,250W,165V Heater harness W Z7524
X3901-81008	AX-65 Dual thermostat harness

Agilent and Competitors' Pumps Cross-Reference Table

Pumping speed m ³ /h (50 Hz)	Agilent Wet Pump	Part Number	Description	Edwards	Pfeiffer	Pfeiffer Adixen	Oerlikon Leybold	Agilent Dry Pump	Part Number	Description
≤ 3	DS 40M	X1699-64052	DS 40M 110V	E2M0.7	DUO 1.6	-	-	IDP-3	IDP3A21	IDP-3 220 V
				E2M1.5	DUO 3					
		X1699-64053	DS 40M 220V	RV3	-	IDP3B21	IDP-3 110 V			
≤ 5	DS 102	9499315M012	DS 102	RV5	DUO 5 DUO 6	Pascal 2005	TRIVAC D4B	IDP-7	X3807-64010	IDP-7
≤ 10	DS 202	9499320M012	DS 202	RV8	DUO 10 DUO 11	Pascal 2010	TRIVAC D8B	IDP-10	X3810-64010	IDP-10
≤ 15	DS 302	9499325M004	DS 302	RV12	-	Pascal 2015	-	IDP-15	X3810-64010	IDP-10
									X3815-64010	IDP-15
≤ 20	DS 402	9499330	DS 402	E1M18 E2M18	DUO 20	Pascal 2021	TRIVAC D16B	IDP-15	X3815-64010	IDP-15
≤ 35	DS 602	9499335M004	DS 602	E2M28	DUO 35	Pascal 2033	TRIVAC D25B	TS 800	X3848-64800	TS 800
40	MS 40	9499225M004	MS 40				SV40Bi			
120	MS 120	X3702-64000	MS 120	-	-	-	-	-	-	-



Contact Your Local Agilent Rep

Phone:
E-Mail:

Call Agilent Toll-Free 00 800 234 234 00
E-Mail: vpt-customer@agilent.com



Order and
track online



Live
Chat

www.agilent.com/chem/vacuum

DE46835373

This information is subject to change without notice.

© Agilent Technologies, Inc. 2023

Published in the USA, October 15, 2023

5994-5449EN